



M COLIBLUE24 Broth ampules (pk/50)

| | |
|-------------------------|----------------|
| Termékszám: | 2608450 |
| HUF Ár (nettó): | Elérhetőségünk |
| Szállítás 2 héten belül | |

In 24 hours, Hach's m-ColiBlue24 Broth give you simultaneous results for total coliforms and *E. coli*. Hach's ready-to-use media eliminates measuring, mixing, and autoclaving steps necessary to prepare dehydrated media. Glass ampules (26084-20), plastic ampules (260845-00), 100mL bottles (26084-42) of m-ColiBlue24 and agar plates (28052-15), afford maximum shelf life and ease of use for pre-measured and prepared media. Enumerate total coliforms and *E. coli* on one petri dish. Differentiate colonies easily - red and blue indicate total coliforms and blue specifies *E. coli*. Get optimal recovery of stressed and injured organisms.

Use m-ColiBlue24 Broth to accurately monitor and evaluate total coliform and *E. coli* in drinking water; wastewater; bottled water; beverages; surface, ground, and well water; and ultrapure, chemical processing and pharmaceutical processing waters.

m-ColiBlue24 should only be used with pads manufactured by the following companies/brands: Pall Corporation (e.g., Sterile Petri Dishes with Absorbant Pads), Private Label/OEM products manufactured by Pall Gelman, Sartorius (cellulose and glass fiber pads), Sartorius BioSart products with the integrated Petri dish, and Millipore EO sterilized pads. No other pad is compatible with the product, and the use of any pad other than those recommended may result in unreliable, inaccurate data.

With ready-to-use m-ColiBlue24 Broth, you can:

Simultaneously detect total coliforms and *E. coli* - in just 24 hours

Simple handling with pre-portioned selective medium

Cost-effective - no fluorescent lamps or colour comparators required

Superior sensitivity down to 1CFU/100 mL

Műszaki adatok

| | |
|------------------------|---------------------|
| Anyag: | Plastic |
| Media Type: | m-Colibblue24 Broth |
| Térfogat/ csomagméret: | 50 /pk |