



Digital titrator cartridge, sodium thiosulfate, 0.200 N

Termékszám: 2267501
HUF Ár (nettó): Elérhetőségünk
Szállítás 2 héten belül

For determination of Dissolved Oxygen by Winkler Method Digital Titrator titration. 0.200 N Na₂S₂O₃. Digits x 0.01 = mg/L DO in 200 mL sample, Digits x 0.00025 = meq.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Méréstartomány:	1 - 10 mg/L DO
Módszer név:	Azide mod. of Winkler Method
Módszer száma:	8215
Paraméter:	Oxygen, Dissolved
Platform:	Digitális titrátor
Tesztek száma:	40 - 100
Titrló:	0.200 N Na ₂ S ₂ O ₃