



Digital titrator cartridge, sodium thiosulfate, stabilized, 0.0250 N

Termékszám:	2409301
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

For dissolved oxygen determination by Winkler titration. Digits x 0.1 = mg/L DO in 25 mL. Stabilized sodium thiosulfate formulation

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Módszer név:	Winkler
Platform:	Oldat