



E.coli and total coliforms, pk/25

Termékszám:	2495525
HUF Ár (nettó):	Elérhetőségünk
Szállítás 2 héten belül	

Presence/Absence Broth with MUG ampules provide an easy-to-use test, ideal for screening drinking water samples for the presence of Total Coliforms and E. coli bacteria. Simply combine the P/A broth with your 100 mL sample and incubate for 24 to 48 hours at 35 °C. A colour change to yellow or yellow brown indicates presence of coliform bacteria.

Complete Solutions for Testing Recreational Water

Complete Solutions for Testing Wastewater

Complete Solutions for Testing Drinking Water

Műszaki adatok

Anyag:	Glass
Térfogat/ csomagméret:	25 /pk