



Digital titrator cartridge, sodium hydroxide, 0.9274 N

Termékszám:

1484201

Veszélyes

ELAVULT CIKK

Ez a termék már nem elérhető.



For volatile acids determination by Digital Titration, 100 - 2400 mg/L Acetic Acid

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Méréstartomány:	100 - 2400 mg/L CH ₃ COOH
Módszer név:	Phenolphthalein End Point
Paraméter:	Volatile Acids
Platform:	Digitális titrátor
Tesztek száma:	40 - 100
Titrló:	0.9274 N NaOH