



Be Right™



Test strips free chlorine, 0-600 mg/L, 100 tests

Termékszám:	2890200
HUF Ár (nettó):	Elérhetőségünk
Elérhető	

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests

Simplify testing and get results in less than 1 minute, in the field or lab

Eliminate chemical handling and clean-up

Reduce guesswork by screening samples

***Not for use in sea or brackish water**

Műszaki adatok

Egység:	100 tests
Lábjegyzet:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Legkisebb beosztás:	Steps: 0, 25, 50, 100, 200, 400, 600
Méréstartomány:	0 - 600 mg/L Cl ₂
Paraméter:	Free Chlorine, High Range
Platform:	Tesz csík