



Silver nitrate digital titrator cartridge, 1.128 N, 1 piece

Termékszám:

1439701

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Elérhető



Ready-to-use titrant solution for the analysis of Chloride (Silver Nitrate Method), 10 - 10000 mg/L Cl⁻.

For chloride determination by Digital Titrator argentometric titration. Hach method 8207. Digits = mg/L Cl⁻ in 50 mL sample. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Méréstartomány:	100 - 10000 mg/L Cl ⁻
Módszer név:	Silver Nitrate
Módszer száma:	8207
Paraméter:	Chloride
Platform:	Digitális titrátor
Tesztek száma:	40 - 100
Titrló:	1.128 N AgNO ₃