



Maintenance kit for 5500 sc Silica Analyser, 2/4 channel

Termékszám:	6788302
HUF Ár (nettó):	Elérhetőségünk
Elérhető	

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Műszaki adatok

A doboz tartalma:	Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions
Product Kit:	5500 sc Silica Analyser, 2/4 Channels

A doboz tartalma

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions