



**Be Right™**



## **Test kit, chlorine, hardness, iron & pH , model CN-39WR**

**Termékszám:**

223002

Veszélyes

HUF Ár (nettó):

## Elérhetőségünk

Szállítási dátum nem látható



Chlorine, Hardness, Iron, and pH Test Kit combines the most popular water treatment test into one compact kit. The chlorine test is USEPA-accepted for drinking water testing if approved by individual states.

**Includes frequently measured water quality parameters in one rugged test kit**

**Uses simple, economical drop count titration and color disc methods**

Tests chlorine (0-3.4 mg/L), hardness (1-20 gpg), iron (0-4 mg/L), and pH (4-10).

## Műszaki adatok

Alkalmazás: Water Quality

## Water Quality

## Méréstartomány:

0 - 3.4 mg/L

## Paraméter

Chlorine H

Platform:

Color

Tesztel

100