



Testkit, glycol visual, model EG-1

Termékszám:

2186400

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Semiquantitative determination of glycol in water or oil. Visual appearance of colour indicates glycol presence.

Ethylene glycol (anti-freeze) contamination of automotive lubricants is an indication of a cracked engine block. Glycol presence also can indicate leaks in cooling systems. Test Kit Model EG-1 presents a purple color if ethylene glycol and other 1,2 glycols are present in oil or water. The kit contains approximately 25 tests and comes with a 14.4" x 8" x 7" Blue Polypropylene Carrying Case.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Méréstartomány:	Presence or absence of glycol, 0- greater than 150 mg/L
Módszer név:	Visual/
	Appearance of color
Paraméter:	Glycol - Ethylene glycol
	in oil or water
Tartó stílus:	G
Tesztek száma:	25
Tesztkészlet model:	EG-1