



Be Right™



## Chlorine (total) test kit, model CN-21P, 10 - 200 mg/L, 100 tests

Termékszám:

2444400

Veszélyes

HUF Ár (nettó):  
Elérhető

Elérhetőségünk



Method / Chemistry: Drop count titration / Thiosulfate

Chlorine Test Kits are manufactured by Hach Company to ensure compliance with stringent quality control standards. Using the Drop Count Titration method, the CN-21P Chlorine Test Kit is used to determine high concentrations of chlorine in sanitizing applications. 10 - 200 mg/L; contains 100 tests.

**Simple, economical titration method**

**Kit contains all required reagents, solutions, and glassware**

### Műszaki adatok

Lábjegyzet:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Méréstartomány:	10 - 200 mg/L
Módszer név:	Drop count titration/
	Thiosulfate
Paraméter:	Chlorine, mid and high range - As total Cl <sub>2</sub>
Platform:	Drop Count Titration
Smallest increment:	10
Tartó stílus:	B
Tesztek száma:	100
Tesztkészlet model:	CN-21P