



## Dissolved oxygen test kit, model OX-2P

**Termékszám:**

**146900**

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható

Veszélyes



Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count titration

Replacement reagent set PN 2439802

**Utilizes drop count titration - a simple, economical titration method**

---

### Műszaki adatok

|                     |                                               |
|---------------------|-----------------------------------------------|
| Méréstartomány:     | 0.2 - 4.0 mg/L & 1 - 20 mg/L O <sub>2</sub>   |
| Módszer név:        | Drop count titration/<br><br>Modified Winkler |
| Paraméter:          | Dissolved Oxygen - As DO                      |
| Platform:           | Drop Count Titration                          |
| Smallest increment: | 0.2                                           |
| Tartó stílus:       | D                                             |
| Tesztek száma:      | 100                                           |
| Tesztkészlet model: | OX-2P                                         |