



Be Right™



Test kit, sulfate, model SF-1

Termékszám:

225100

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



Widely found in natural waters, sulfate may be at high levels in mine drainage. High concentrations of magnesium or sodium sulfate in water act as cathartics. This test kit provides economical determination of sulfate levels using the Extinction/Turbidimetric methods.

Economical test kit for determination of sulfate levels using the Extinction/Turbidimetric Method

Műszaki adatok

Méréstartomány:	50 - 200 mg/L
Módszer név:	Extinction/Turbidimetric
Paraméter:	Sulfate - As SO ₄ ²⁻
Smallest increment:	50
Tartó stílus:	D
Tesztek száma:	100
Tesztkészlet model:	SF-1