



Test kit, alkalinity, model AL-AP

Termékszám:

2444300

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Lábjegyzet:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. **1 gpg = 17.1 mg/L (MO) = Methyl Orange
Méréstartomány:	0 - 8 gpg & 1 - 20 gpg
Módszer név:	Drop count titration/Sulfuric acid
Paraméter:	Alkalinity - (P) & Total (MO) As CaCO ₃
Platform:	Drop Count Titration
Smallest increment:	0.4 gpg
	1 gpg
Tartó stílus:	B
Tesztek száma:	100
Tesztkészlet model:	AL-AP