



Be Right™



## Alkalinity test kit AL-AP, 5 - 100 / 20 - 400 mg/L CaCO<sub>3</sub>, 100 pcs

Termékszám:

2444301

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists.

Method: Drop count titration

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Kit utilizes drop count titration method and contains 100 tests.

### Simple, economical titration method

Rugged kit contains all required reagents, solutions, and glassware

### Műszaki adatok

Lábjegyzet:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Méréstartomány:	5 - 100 mg/L & 20 - 400 mg/L
Módszer név:	Drop count titration/Sulfuric acid
Paraméter:	Alkalinity - (P) & Total (MO) As CaCO <sub>3</sub>
Platform:	Drop Count Titration
Smallest increment:	5
	20
Tartó stílus:	B
Tesztek száma:	100
Tesztkészlet model:	AL-AP, MG-L