



Be Right™



Stablcal Zavarosság standard kalibráló készlet, 0 - 4000 NTU

Termékszám:

2662105

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Elérhető



5 használatra kész, stabilizált formazin elsődleges standard oldatok. Készlet légmentes fiolákat tartalmaz az alábbiakhoz : <0.1 NTU, 20 NTU, 200 NTU, 1000 NTU, 4000 NTU.

Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

Stablcal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilised version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

Műszaki adatok

Koncentráció: < 0,1, 20, 200, 1000, 4000 NTU

Műszer: 2100N, 2100N IS, TL2300 and TL2310 IS

Paraméter: Turbidity

Térfogat/ csomagmáret: Set of 5