



Reagent set, hardness (total), digital titrator, full range

Termékszám:

2272000

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



For the determination of Total Hardness, 10-160 mg/l CaCO₃. EDTA Method.

For Total Hardness determination by EDTA Digital Titrator titration with Calmagite indicator. Range: 10-160 and 100-4000 mg/L as CaCO₃. Set contains ManVer® 2 Indicator, Hardness 1 Buffer Solution and one cartridge each of 0.0800 M and 0.800 M EDTA.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Műszaki adatok

Méréstartomány:

10 - 160 mg/L & 100-4000 mg/L CaCO₃

Módszer név:

EDTA Titration

Paraméter:

Hardness (Total)

Speciális tulajdonság:

Reagent set, Digital Titrator