



Nátrium-klorid, 10 000 µS/cm, 500 ml

Termékszám:	2972249
HUF Ár (nettó):	Elérhetőségünk
Szállítás 2 héten belül	

Conductivity standard solutions for laboratory and field analyses

Conductivity standard solution, 10000 µS/cm $\pm$ 100 µS/cm at 25°C, 5.556 g/L $\pm$ 0.03 g/L as Sodium Chloride (NaCl), (1 µS/cm $\sim$ 0.5 mg/L TDS), 500 mL bottle.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Műszaki adatok

Koncentráció:	10000 µS/cm at 25°C $\pm$ 100 µS/cm at 25°C
Leírás:	Nátrium-klorid, 10 000 µS/cm, 500 ml
Package settings:	Bottle
Paraméter:	Vezetőképesség
Platform:	Oldat
Tárolási feltételek:	10 - 25 °C (protect from light)
Térfogat/ csomagméret:	500 mL