



HQ11D Digital pH meter kit, pH Gel electrode, Outdoor, 5m

Termékszám:

HQ11D.99.105000

ELAVULT CIKK

Ez a termék már nem elérhető.

Digital Electrochemistry with HQD

Portable pH/ORP meter. Instrument with Outdoor pH Gel electrode, cable length 5m.

The digital meter/electrode system combines reliability, flexibility and ease of use. Interchangeable Intellical electrodes are recognised automatically and store all relevant data. Rugged, nearly indestructible outdoor versions with multiple cable lengths enable measurements even at previously inaccessible locations.

Replaced by NEW HQ Series.

Outstanding reliability and simple handling

Versatile electrodes for all applications, e.g. wastewater, drinking water, process water

Reliable results from inaccessible measurement locations and over long distances - even for pH

Full GLP data management

Műszaki adatok

A doboz tartalma:	Meter package includes HQ11D Portable Meter; 4 AA batteries; PHC101 Field Rugged Gel pH Electrode, 5 m cable; Singlet pH 4 & 7 pH Buffer 10 of each; quick-start guide and user manuals.
Adat export:	Download via USB connection to PC or flash memory device. Automatically transfer entire data log or as readings are taken. Power supply adapter required.
Adattárolás:	Automatic, GLP ISO compliant reading data stored with batteriescalibration details.
Adattárolás/memória:	500 records/FIFO
Áram és USB működés:	Optional
Auto-puffer felismerés:	IUPAC standards (DIN 19266), Technical buffer (DIN 19267), 4-7-10 series, or customer defined
Bemenetek:	1
Felhasználói felület nyelve:	English, German, French, Italian, Spanish, Danish, Dutch, Polish, Portuguese, Turkish, Finnish, Czech, Russian.
Fix puffer választás:	IUPAC standards (DIN 19266), Technical buffer (DIN 19267), 4-7-10 series, or customer defined
Garancia:	36 hónap
GLP jellemzők:	Date; Time; Sample ID; Operator ID
Hőmérséklet felbontás:	0,1 °C

Hőmérséklet kompenzáció:	A pH automatikus hőmérséklet-kompenzációja
Hőmérséklet mérés:	-10 - 110 °C
Hőmérséklet pontosság:	± 0.3 °C
Includes:	Meter package includes HQ11D Portable Meter; 4 AA batteries; PHC101 Field Rugged Gel pH Electrode, 5 m cable; Singlet™ pH 4 & 7 pH Buffer 10 of each; quick-start guide and user manuals.
Járna-e hozzá szondák?:	PHC10105
Kalibrációs intervallum/emlékeztető/riasztás:	Ki, 2 órától 7 napig választható
Készlet?:	Igen
Kijelző:	Detailed mode/Large mode
Kijelző típus:	240 x 160 pixel LCD with batteriesbacklight illumination
Kijelző zár funkció:	Continuous / Auto-stabilization ("press to read") / At Interval
Környezeti feltétel: hőmérséklet:	0 - 60 °C (32 - 140 °F)
Környezeti feltétel: relatív páratartalom:	90 % (nem kondenzáló)
Megfelelőségi igazolások:	CE.WEEE
Mérés módszer:	Probe specific programmed method settings
Működési hibaüzenet:	Clear text error messages displayed
Működési interfész:	Soft Touch Keypad
Műszer:	Portable
mV felbontás:	0.1 mV
mV mérés :	-1500 - 1500 mV
mV mérés stabil leolvasásnál:	Yes
Paraméter:	pH/Oxydo Reduction Potential (ORP)
Parameters:	pH
	NA
	NA
PC data transfer software :	HQD Series Meter Data Transfer Utility
pH Electrode calibration:	1 - 3 Calibration points
	Calibration summary data logged and displayed
pH felbontás:	Választható funkciók 0.001 - 0.1 pH
pH Measurement Range:	0 - 14 pH
Printer:	Yes, Optional
Redox elektróda kalibráció:	Predefined ORP standards (including Zobell's solution)
Sensor A:	PHC10105 (5 m rugged)
Sensor B:	NA
Sensor C:	NA
Sensors:	Sensor A: PHC10105
	Sensor B: NA
	Sensor C:
Súly:	338 g (0.75 lb) without batteries; 430 g (0.95 lb) with batteries
Szokásos kalibrációs standard:	Yes
Szondatípus:	Intellical Digital SMART Standard Laboratory or Rugged Field
Tartalom:	No

Tartó IP besorolás:	IP67
Típus:	HQ11D – pH/1 Channel
Type of measurement:	Field Measurements

A doboz tartalma

Meter package includes HQ11D Portable Meter; 4 AA batteries; PHC101 Field Rugged Gel pH Electrode, 5 m cable; Singlet pH 4 & 7 pH Buffer 10 of each; quick-start guide and user manuals.