



DREL Portable Industrial Water Lab DR1900

Termékszám:

LZV966

ELAVULT CIKK

Ez a termék már nem elérhető.

Portable water lab including spectrophotometer DR1900, complete with batteries and cuvettes, instrument case, reagent & apparatus case, reagent set, and instrument set.

Today, we have to be able to carry out analysis anywhere. The portable laboratory contains everything needed to ensure that accurate measurements can be taken.

Even in the field, there is the need for efficient measurement technology, practical accessories and ready-to-use reagents. The integrated work surface enables analysis to be carried out safely and easily at any location.

Complete water lab for field testing, incl. photometer, digital titrator and electrochemical methods

Versatile and rugged

Simple operation with ready-to-use reagents

Műszaki adatok

Garancia:	24 hónap
Hullámhossz tartomány:	340 nm - 800 nm
Includes:	DR 1900 portable Spectrophotometer (DR1900-02L), instrument case (5838500), reagent case (5839400), reagent set (2923300), reagent set (2923400)
Méreték (MxSZxM):	280 mm x 533 mm x 290 mm