



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Termékszám:	108076
HUF Ár (nettó):	Elérhetőségünk
Szállítás 2 héten belül	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Műszaki adatok

Anyag:	Polypropylene, Low Form
Anyaga:	Polypropylene
Átmérő, külső:	83.6 mm with spout
Egység:	6 /pk
Kapacitás:	250 mL