



Nitrate Standard Solution, 1000 mg/L, 500 mL

Termékszám:

1279249

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Elérhető



Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Nitrate, 1000 mg/L ± 10 mg/L $\text{NO}_3\text{-N}$ (NIST). 500 mL bottle

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Műszaki adatok

Koncentráció:	4430 $\pm$ 44 mg/L NO_3
Különleges jegyzet:	Observed shelf life after opening, following recommended storage and use: 3 months
Leírás:	Nitrogen-Nitrate Standard Solution as $\text{NO}_3\text{-N}$ (NIST)
Package settings:	Bottle
Paraméter:	Nitrate
Platform:	Oldat
Tárolási feltételek:	10 - 25 °C (protect from light)
Térfogat/ csomagméret:	500 mL