



Buret, Acrylic, 50 mL

Termékszám:	1405941
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Durable, crystal-clear acrylic body is individually calibrated to meet the accuracy requirements of ASTM E287. Leak-proof stopcock with self-lubricating PTFE TFE plug. Excellent for general lab or field use in schools, industrial, municipal or food applications. Unaffected by dilute mineral acids and bases (except hydrofluoric acid and ammonium hydroxide). Can contain liquids within a pH range of 1 - 14 and up to 5% wt/wt of acids or 30% bases can be used. Not designed for exposure to alcohol or organic solvents.

Műszaki adatok

Anyaga:	Acrylic
Graduations:	0.1 mL
Kapacitás:	50 mL
Speciális tulajdonság:	Buret