



## Standard solution, alkalinity, 25000 mg/L as CaCO<sub>3</sub> (NIST), pk/16 - 10 mL Voluette ampules

|                    |                |
|--------------------|----------------|
| <b>Termékszám:</b> | <b>1427810</b> |
| HUF Ár (nettó):    | Elérhetőségünk |
| Elérhető           |                |

For accuracy checks in alkalinity determinations

### Simple accuracy checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the Tensette pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

### Minimise waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

### Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

---

## Műszaki adatok

|                        |                                       |
|------------------------|---------------------------------------|
| Koncentráció:          | 0,5 N (25000 mg/L CaCO <sub>3</sub> ) |
| Paraméter:             | Alkalinity                            |
| Platform:              | Oldat                                 |
| Térfogat/ csomagméret: | 16 x 10 mL ampules                    |