



**Be Right™**



## **Manganese test kit, model MN-5**

**Termékszám:**

146700

## Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. An economical solution, this kit uses the simple, accurate color disc method.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

## Műszaki adatok

Méréstartomány:	0 - 3.0 mg/L Mn
Módszer név:	Color disc/Cold periodate
Paraméter:	Manganese - As Mn
Parameter:	Manganese
Platform:	Color Disc
Smallest increment:	0.1
Standard mérési tartomány:	0 - 3.0 mg/L Mn
Tartó stílus:	D
Tesztek száma:	100
Tesztkészlet model:	MN-5
Típus:	MN-5