



Be Right™



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 500 mL, 12/pk

Termékszám:	1472459
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Műszaki adatok

Anyag:	PP
Body Material:	Plastic
Kapacitás:	500 mL
Package size UOM:	12pk
Térfogat/ csomagméret:	12 /pk