



**Be Right™**



## Hydrazine color disc test kit, model HY-2

**Termékszám:**

**184900**

**Veszélyes**

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method with a range of 0.02-1 mg/L, this kit contains 300 tests.

**Measures Hydrazine, a common oxygen scavenger in boiler water, using the simple color disc method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Műszaki adatok

Méréstartomány:	0 - 1 mg/L N <sub>2</sub> H <sub>4</sub>
Módszer név:	Color disc/  p-Dimethylamino-  benzaldehyde
Paraméter:	Hydrazine - As N <sub>2</sub> H <sub>4</sub>
Parameter:	Detergents
Platform:	Color Disc
Smallest increment:	0
Standard mérési tartomány:	0 - 1 mg/L N <sub>2</sub> H <sub>4</sub>
Tartó stílus:	D
Tesztek száma:	300
Tesztkészlet model:	HY-2
Típus:	HY-2