



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 500 mL

Termékszám: 2089649

ELAVULT CIKK

Ez a termék már nem elérhető.

Javasolt cserék:

- 2089679

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Műszaki adatok

Anyag:	Glass Screw Cap
Egység:	1 each
Kapacitás:	500 mL