



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 250 mL

Termékszám:

2089676

ELAVULT CIKK

Ez a termék már nem elérhető.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Műszaki adatok

Anyag:	Glass Screw Cap
Egység:	6/pk
Kapacitás:	250 mL