



## Reagent and apparatus set, water in oil

**Termékszám:**

**2237300**

**Veszélyes**

**ELAVULT CIKK**

**Ez a termék már nem elérhető.**



### Model WO-1 Test Kit

The Model WO-1 test kit provides critical data in oil quality determinations, oil line operation, and diesel engine maintenance. The kit has proven particularly useful on board ships. Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated.

**Measures low levels of water in oil**

**Provides critical data in oil quality deterinations, oil line operation, and diesel engine maintenance**

**Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated**

---

### Műszaki adatok

Méréstartomány:	0.05 - 1 %
Méréstartomány (2):	0 - 10 %
Módszer név:	Volumetric displacement/ Calcium hydride
Paraméter:	Water-in-oil - As %
Smallest increment:	0.05% 0.20%
Tartó stílus:	no case
Tesztek száma:	25
Tesztkészlet model:	WO-1