



SA filter, cone, diameter 110 mm, soil analysis grade, 500/pk

Termékszám:	2263356
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Designed for sandy soils, SA filters are conically folded and have been optimized for use with the common extraction solutions used in soil nutrient analysis (Mehlich 1, Mehlich 3, Morgan, Modified Morgan, Bray, and Olsen). SA Filters provide the best combination of retention rating, filtration speed and purity requirements for soil nutrient testing (P, K, Ca, Mg). Medium filter retention.

Műszaki adatok

Anyag:	Cellulose
Átmérő:	110 mm
Térfogat/ csomagméret:	50 /pk