



## A-1 broth medium, dehydrated, 500 g

**Termékszám:**

**2309934**

**Veszélyes**

**ELAVULT CIKK**

Ez a termék már nem elérhető.



A-1 Broth Medium, Dehydrated. Broth for 24-hour MPN method for fecal coliforms. See Microbiological Testing for prepared medium. Bottle/500 g.

### **Complete Solutions for Testing Recreational Water**

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

### **Complete Solutions for Testing Wastewater**

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

### **Complete Solutions for Testing Drinking Water**

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

---

## **Műszaki adatok**

Media Type: A-1 Broth Medium

Térfogat/ csomagméret: 500 g