



Filter holder, for vacuum filtration, glass

Termékszám:	234000
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Beaker-style funnel (300 mL) with 47 mm membrane support for filtration of aqueous samples. Fritted borosilicate glass membrane support for even filtration across surface of membrane. Funnel attaches evenly to base via aluminum clamp. Filter flask sold separately.

Műszaki adatok

A doboz tartalma:	Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).
Átmérő:	47 mm
Térfogat/ csomagméret:	1

A doboz tartalma

Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).