



Be Right™



Test kit, oxygen scavenger, color disc, model OS-1, 50 tests

Termékszám:

2349500

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Műszaki adatok

Méréstartomány:	0 - 52 µg/L DEHA (and other Oxygen Scavengers)
Méréstartomány (2):	0 - 1360 µg/L DEHA (and other Oxygen Scavengers)
Módszer név:	Color disc/Iron reduction
Paraméter:	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Platform:	Color Disc
Tartó stílus:	R
Tesztek száma:	50
Tesztkészlet model:	OS-1
Típus:	OS-1