



Filter, membrane, black, pore size 0.45 µm, diameter 47 mm, 200/pk, w/out pads

Termékszám:	2490300
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Provides contrast for microbiological and particulate analysis. Modified polyethersulfone creates dark background for excellent contrast useful in enumerating opaque colonies and light-covered particulate. White grid line aids enumeration. Certified for membrane filtration analysis. Sterile and individually packaged. Pore size 0.45 µm. Diameter: 47 mm.

Műszaki adatok

Anyag:	Modified polyethersulfone
Átmérő:	47 mm
Filter Type:	Filter
Szemcseméret:	0.45 µm
Térfogat/ csomagméret:	200 /pk