



Be Right™



Bottle, sampling, presterilized, round, 12/pk

Termékszám:

2495012

ELAVULT CIKK

Ez a termék már nem elérhető.

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Works well with P-A Broth Ampules and PathoScreen Medium.

Műszaki adatok

Dechlorinating Agent?:

No

Térfogat/ csomagméret:

12/pk