



Test Kit, AccuVac® Color Disc, Ozone, 0.05-1.50 mg/L

Termékszám:

2518050

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



This ozone test kit uses the fast, accurate visual method for determining ozone concentrations. Use to control disinfectant levels in water, wastewater, and industrial water. Kit contains 25 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

AccuVac® vials mean no mess, no measuring

Indigo trisulfonate method masks interferences from chlorine and other oxidizing agents and makes ozone analysis easier

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Műszaki adatok

Méréstartomány:	0 - 1.50 mg/L O ₃
Módszer név:	Color disc/ Indigo trisulfonate
Paraméter:	Ozone - As O ₃
Parameter:	Ozone
Platform:	Accuvac
Smallest increment:	0.05
Standard mérési tartomány:	0 - 1.50 mg/L O ₃
Tartó stílus:	G
Tesztek száma:	25
Tesztkészlet model:	AccuVac® Kit
Típus:	NA