



Filter, TSS glass fiber, pore size 1.5 µm, diameter 110 mm, 100/pk

Termékszám:	2569100
HUF Ár (nettó):	Elérhetőségünk
Szállítás 2 héten belül	

TSS Filters are recommended for use in Standard Methods 2540 C and 2540 D. Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results. Fast flow rate, fine porosity and 1.5 µm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater.

Műszaki adatok

Anyag:	Glass Microfiber
Átmérő:	110 mm
Szemcseméret:	1.5 µm
Térfogat/ csomagméret:	100 /pk