



## Rack, test tube, no-wire

|                              |                |
|------------------------------|----------------|
| <b>Termékszám:</b>           | <b>2587300</b> |
| HUF Ár (nettó):              | Elérhetőségünk |
| Szállítási dátum nem látható |                |

Ninety-place No-Wire Rack features locking mechanism that prevents “floaters” in water baths and speeds decanting by allowing the whole rack to be inverted. Heavier-than-water polypropylene assures non-floating racks when immersed, even in agitated water. Autoclavable. Numbered and lettered rows (6 x 15 array) aid in identification of tubes (10-13 mm diameter). Collapsible for storage, yet assemble in seconds. Overall dimensions: 24.8 x 10.5 x 7 cm (9.75 x 4.12 x 2.75").

---

### Műszaki adatok

|                        |             |
|------------------------|-------------|
| Hossz:                 | 248 mm      |
| Mélység:               | 2.75 Inches |
| Szélesség:             | 4.12 Inches |
| Térfogat/ csomagméret: | 1 each      |