



Test kit, RapidSilver[™] visual

Termékszám: 2674500

ELAVULT CIKK
Ez a termék már nem elérhető.

Test Kit Silver, Model RapidSilver[™]. Range 0-50 ppb. Silver As Ag Visual Comparison.

Provides fast, accurate information on silver discharges

Used by surface finishers, photographic film manufacturers and processors, and others to prevent pollution and economic loss

Kit directions also include digestion and dilution procedures needed for testing silver strike bath samples

Műszaki adatok

Méréstartomány:	5 - 50 ppb Ag
Módszer név:	Visual Comparison
Paraméter:	Silver-As Ag
Parameter:	Silver
Smallest increment:	5 µg/L
Standard mérési tartomány:	5 - 50 ppb Ag
Tartó stílus:	F
Tesztek száma:	100
Tesztkészlet model:	RapidSilver [™]
Típus:	NA