



CO2 calibration starter kit, 0-3g/kg

Termékszám: 32315

ELAVULT CIKK
Ez a termék már nem elérhető.

CO2 calibration starter kit, 0-3g/kg. Includes calibration bottle, adapter, and reagents from part number 32316. Suitable for ORBISPHERE instruments employing Thermal Conductivity sensors.

Traceable and Trustworthy Results

Effective calibration is critical to the performance of any analytical instrument and the ORBISPHERE 3231X is designed to give a traceable and trustworthy calibration standard with defined carbon dioxide content to +/- 0.5 %. Each reagent kit contains a certificate of analysis.

Quick and Easy to Use

Műszaki adatok

Méréstartomány:	0 - 3 g/kg
Műszer:	Orbisphere 3658 Analyzer Portable CO2
Paraméter:	CO2
Tárolási feltételek:	Store in a dry place at room temperature