



## Large custom kit case assembly (574 x 325 x 295 mm)

|                              |                |
|------------------------------|----------------|
| <b>Termékszám:</b>           | <b>4766600</b> |
| HUF Ár (nettó):              | Elérhetőségünk |
| Szállítási dátum nem látható |                |

Designed to withstand the harshest conditions, this rugged, blue polypropylene case assembly includes a foam insert to prevent kit components from shifting, breaking or spilling.

**Allows easy transport of instruments, reagents, and apparatus into the field**

**Chemical-resistant test kit case is molded from durable, easy-to-clean polypropylene**

---

### Műszaki adatok

|                   |                          |
|-------------------|--------------------------|
| Méretek (MxSZxM): | 574 mm x 325 mm x 295 mm |
|-------------------|--------------------------|