



Free Ammonia Reagent Set, 50 tests

Termékszám:

2879200

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Elérhető



Improved reagent set for sequential determination of free ammonia and monochloramine by the Indophenol method. **Replaces 2879700.**

Requires new procedure (See Reagent/Procedure tab) and instrument software update (See Downloads>Other>Product Notice-Monochloramine and Free Ammonia).

Hach Method 10200 Ranges: 0.04 - 4.50 mg/L Cl_2 and 0.02 - 0.50 mg/L $\text{NH}_3\text{-N}$. Set contains reagents for approximately 50 tests.

Műszaki adatok

Méréstartomány:	0.02 - 0.50 mg/L $\text{NH}_3\text{-N}$
Módszer név:	Indophenol
Módszer száma:	10200
Műszer:	All except DR/820
Tesztek száma:	50