



Xenomax TSS-FVS Filter 47 mm, 100/pk

Termékszám:	2511307
HUF Ár (nettó):	Elérhetőségünk
Szállítás 2 héten belül	

Xenomax TSS-FVS Filters are specifically designed to ensure compliance with the filter preparation requirements of Standard Methods 2540 D and Standard Methods 2540 E (Fixed and Volatile Solids Ignited at 550 °C). Independent 3rd party evaluation of Xenomax and 934-AH filters in both methods give demonstrably equivalent results.

In addition, these 1.5 µm binder-free, fast flowing, fine porosity borosilicate glass microfiber filters have been optimized for better wet strength/handleability and high temperature (550 °C) blank stability. FVS Filters are available in popular diameter sizes and ready for use in a wide variety of labware including Buchner funnels and Gooch Crucibles.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare.

Műszaki adatok

Anyag:	Glass Microfiber
Átmérő:	47 mm
Szemcseméret:	1.5 µm
Térfogat/ csomagméret:	100 /pk