



**Be Right™**



## EZ3005 Chloride Analyser

**Termékszám:**

**EZ3005.XXXXXXXX**

**HUF Ár (nettó):**

**Elérhetőségünk**

### Online ion-selective Chloride measurements in industrial and environmental applications

#### ISE technology for optimal analytical performance

With limited maintenance requirements and reduced reagent consumption, the EZ3000 Series are the ideal choice for a wide range of water monitoring applications where ion-selective electrodes are the preferred analytical technique. Outstanding precision and stability is guaranteed by the temperaturecontrolled measurement.

#### Direct, discontinuous ISE method

Contrary to separate electrodes or other analysers in the market, the EZ3000 Series does not run a continuous measurement. The principle of discontinuous ISE analysis not only enhances control over conversion of ion activity to electric potential, it also eliminates risk of cross-contamination between cycles and reduces overall consumption of reagents.

The EZ3000 Series combine unique ISE technology with a set of analysis, control and communication features in an industrial analyser mainframe with designed for the highest performance:

- Automatic direct ion-selective measurements
- Smart automatic features
- Control and communication via industrial panel PC
- Standard 4 - 20 mA signal output with alarm processing
- Communication ports supporting connectivity to Modbus
- Multiple stream analysis
- Reduced reagent consumption

There are many additional options available. Please contact Hach for more details.

#### Műszaki adatok

A doboz tartalma:

EZ3005 Chloride Analyser, Instruction Manual, 1 x Double Bit Door Key, 1 x Mounting Brackets, 1 x ISE Electrode and Chloride Electrode, combined, and 1 x empty 2.5L Reagent Container with Fittings (for Buffer Solution)

Analóg kimenetek:

Active 4 - 20 mA max. 500 Ohm load, standard 1, max. 8 (option)

Anyag:	Hinged part: Thermoform ABS, door: PMMA
Automatikus tisztítás:	Wall section: Galvanised steel, powder coated
Demineralised water:	Igen
Detektálási határ:	For rinsing
Digitális kimenetek:	25 mg/L
Earth connection:	Optional: Modbus (TCP/IP, RS485)
Garancia:	Dry and clean earth pole with low impedance (< 1 Ohm) using an earth cable of > 2.5 mm <sup>2</sup>
Instrument air:	24 hónap
Interferenciák:	Dry and oil free according to ISA-S7.0.01-1996 quality standard for instrument air
Kalibrálás:	Bromide, sulphide, iodide, cyanide ions may interfere. Mercury must be absent. Ammonia and thiosulphate may interfere. Fats, oil, proteins, surfactants and tar.
Környezeti hőmérséklet:	Automatic, 2-point; frequency freely programmable
Lefolyócső:	10 - 30 °C ±4 °C deviation at 5 - 95% relative humidity (non-condensing)
Mérés módszer:	Atmospheric pressure, vented, min. Ø 32 mm
Méréstartomány:	Discontinuous, direct measurement by combined ion-selective electrode, conform with standard methods EPA 9212 and ASTM D512-12
	100 - 1,000 mg/L Cl <sup>-</sup>
Méretek (MxSzxm):	Optional:
Minta hőmérséklet:	25 - 250 mg/L Cl <sup>-</sup>
Minta minőség:	690 mm x 465 mm x 330 mm
Minta: nyomás:	10 - 30 °C
Mintaáram:	Maximum particle size 100 µm, < 0.1 g/L; Turbidity < 50 NTU
Mintacellák:	By external overflow vessel
Mintautak száma:	100 - 300 mL/min
Módszer név:	1 minutes
Paraméter:	2, 4 or 8
Pontosság:	Ionselective
Reagent Requirements:	Chloride
Riasztás:	Better than 2% full scale range for standard test solutions
Súly:	Keep between 10 - 30 °C
Szünet mód:	1x malfunctioning, 4x user-configurable, max. 24 VDC/0.5 A, potential free contacts
Tanúsítványok:	25 kg
Teljesítmény:	CE compliant / UL certified
Típus:	100 - 240 VAC, 50/60 Hz
Validation:	100 - 240 VAC, 50/60 Hz
Védeeltségi osztály:	Max. power consumption: 120 VA
	EZ3005
	Automatic; frequency freely programmable
	Protection class: Analyser cabinet: IP44 / Panel PC: IP65

## A doboz tartalma

EZ3005 Chloride Analyser, Instruction Manual, 1 x Double Bit Door Key, 1 x Mounting Brackets, 1 x ISE Electrode and Chloride Electrode, combined, and 1 x empty 2.5L Reagent Container with Fittings (for Buffer Solution)