



Be Right™



Test strips total hardness, 0-425 mg/L

Termékszám: 2745250
HUF Ár (nettó): Elérhetőségünk
Elérhető

Test strips in a variety of parameters.

Water hardness is caused almost entirely by calcium and magnesium ions. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. Contains 50 tests.

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Műszaki adatok

Lábjegyzet:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Legkisebb beosztás:	Steps: 0, 25, 50, 120, 250, 425 ppm
Méréstartomány:	0 - 425 ppm
Méréstartomány (2):	0 - 25 gpg
Paraméter:	Hardness, Total As CaCO ₃
Platform:	Tesz csík
Smallest increment:	Intervalos de: 0, 1.5, 3, 7, 15, 25 gpg
Szállítási tömeg:	0.2
Tesztek száma:	50
Tesztkészlet model:	Test Strips