



Be Right™



## Funnel, Buchner, Polypropylene, 186 mL Approximate Volume

**Termékszám:** 2090553  
**HUF Ár (nettó):** Elérhetőségünk  
**Szállítás 2 héten belül**

Light weight and durable polypropylene construction with a removable filter plate for easy cleaning, chemical resistant. Fully autoclavable. Stopper sold separately.

Capacity: 186 mL

Inside Top Diameter: 76 mm

Outer Diameter of stem: 11 mm

Stem Length: 71 mm

Overall Height: 146 mm

Fits Filter Paper Size: 70 mm

Autoclavable: up to 135 °C

---

### Műszaki adatok

Átmérő:	11 mm (Stem)
Hossz:	Stem: 65 mm
I. D.:	76 mm
Kapacitás:	186 mL
Térfogat/ csomagméret:	1 each