



## Soxhlet thimbles, cellulose, 33 x 80 mm, 25 pk

**Termékszám:**

**2966700**

**ELAVULT CIKK**

Ez a termék már nem elérhető.

Xenosep extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

Xenosep extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

The thimbles are manufactured from high purity alpha cellulose cotton linters which provide excellent mechanical strength, stability, and retention with very low/no background extractables in a disposable, sample holder format.

---

### Műszaki adatok

Anyag:	Cellulose
Szemcseméret:	10 µm
Térfogat/ csomagméret:	25 /pk