



Colony Counter, Quebec, 220 Vac, 50 Hz

Termékszám: 2252102

ELAVULT CIKK

Ez a termék már nem elérhető.

For counting bacterial colonies. Features 11.5 cm (4.5") lens with 1.5x magnification to eliminate parallax errors; 40-watt tungsten bulbs to illuminate dark backgrounds; Wolffhuegel guide plate and 3-wire cord. Accommodates Stewart plates. Shee metal case. Dimensions: 25 x 28 x 27 cm (10 x 11 x 10.5"). CSA certified.

Illumination

Illumination is provided by 40W tungsten bulbs.

Magnification

A 113mm lens offers 1.5x magnification and is positioned to eliminate parallax errors.

Műszaki adatok

A doboz tartalma: Includes: Colony counter, guide plate, Instructional manual and power cord.

Áramellátási követelmények (feszültség): 220 V

A doboz tartalma

Includes: Colony counter, guide plate, Instructional manual and power cord.