



Iron Reagent Powder Pills, 0.02-3.00mg/L Fe, FerroVer

Termékszám:

85428

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Szállítás 2 héten belül



Powder reagent in sealed package, for 25mL sample.

FerroVer Iron Reagent: For the determination of total soluble Iron including all ferric (Fe_3^+) and ferrous (Fe_2^+) Iron as well as complexed Iron compounds (Fe-EDTA). Perfect for many applications and a wide range of concentrations. Determines Iron in the presence of Copper.

Műszaki adatok

EPA compliant:	Igen
Méréstartomány:	0,02 - 3,00 mg/L Fe
Módszer név:	FerroVer
Módszer száma:	8008
Platform:	Porpárna
Térfogat/ csomagméret:	1000 /pk
Tesztek száma:	1000